

PATIENT BENEFITS

- Affordable
- Comfortable
- Painless
- Non-invasive application
- Wide treatment range
- Easily replaced
- Good dental hygiene
- Interactive



DESIGNED WITH DENTISTS
FOR DENTISTS

SCAN ME FOR



DIGITAL BROCHURE



Medicore
CLEAR ORTHODONTIC ALIGNERS

DENTIST BENEFITS

- 10 Day turn around time
- Manufactured in RSA
- Digital marketing for your practice
- Interactive online platform
- Designated sales consultant



HOW THE PROCESS WORKS

CASE ENTRY:

Step 1 – Panoramic X-ray

Easily upload a Panoramic X-ray using our Medicore Platform.

Step 2 – Clinical prescription form

Medicore's compact Clinical Prescription form is used to convey special instructions the Orthodontist.

Step 3 – Clinical photographs:

Easily upload your clinical photos using our Medicore platform.

Step 4 – Digital Scan:

Upload the patients scan directly from your scanner. Our team offers on-site training and installation.

MEDICORE PROCESS:

Step 1 – Digitalisation:

A digital simulated treatment plan is sent to the dentist for approval before manufacturing starts.

Step 2 – Manufacturing:

Manufacturing of Medicore Clear Aligners is done using high-tech 3D-Printing technology and registered dental technicians.


Step 3 – Delivery:

Delivery is available country-wide!

UNMATCHED QUALITY

Material- CA Pro

High quality material ensures thinner less invasive aligners. The three layer aligner material used maintains its rigidity throughout the treatment period.



Tolerance – 25 micron 3D printing

Medicore uses state of the art 3D printers that produce dental models 30% more accurate than conventional methods.

